Tech Sheet - Tech Device Tattoos

Size Available: 1 standard size with 8 available formats, custom sizes and formats available

Process: Full Color Digital Print

1-sided printing available

Vertical or Horizontal Orientation Kiss-cut or Straight Cut available

Material: Full Color Gloss Tech Device Tattoos

4 mil white, polymeric calendared vinylGrey, removable solvent-based adhesive

- 7 mil. Double PE-coated silicone paper, air-release liner designed for bubble-free installation
- · Gloss, smooth finish
- Removable and Repositionable
- · For use on smooth, flat and slightly curved surfaces
- · Indoor use only
- Minimum temperature for application is 40 degrees Fahrenheit
- Anticipated life span of product is 3-5 years for indoor applications.

Full Color Matte Tech Device Tattoos - Not a Standard Offer

- · 6 mil white, polymeric calendared vinyl
- · Microsphere acrylic adhesive
- 90# layflat liner
- · Matte, smooth finish
- Removable and Repositionable
- · For use on smooth, flat and slightly curved surfaces
- · Indoor use only
- Minimum temperature for application is 40 degrees Fahrenheit
- Anticipated life span of product is 1+ year (s) for indoor applications.

Full Color Fabric Tech Device Tattoos

- 6 mil. white, polyester fabric film made from 100% post-consumer recycled fabric, made in the USA
- Acrylic adhesive
- · 2 mil. Peel-able paper liner
- Matte, slight texture finish
- · High opacity material that eliminates show-through
- Removable and Repositionable
- For use on smooth, flat and slightly curved surfaces
- · Indoor/Outdoor use
- Minimum temperature for application is 40 degrees Fahrenheit
- Anticipated life span of product is 5+ year (s) for indoor applications. Life span significantly decreases to 6 months to 1 year depending on environment.
- · Not recommended for brick, non-painted cement surfaces, or heavy painted spackled wall board.
- Environmentally safe and biodegradable over time, no PVC, is non-toxic.

Production Time: 2 working days

*Quantities of 50 or more and custom sizes may require additional production time

Set-up Charges: FREE

Art Charges: FREE Mino

Minor artwork adjustments can be made with no additional fees applied. This would include items such as slight modification to size and/or color, name/text changes, and artwork clean up (i.e. removing spots from scanned artwork). These alterations should take no longer than 10-15 minutes to complete. Additional design/layout/color assistance is available for an hourly rate. This would include items such as color matching or color separation, logo redraw/manipulation or any additional layout design work needed on artwork.

Standard Offering Sizes:


Gloss Tech Device Tattoo Item Number	Fabric Tech Device Tattoo Item Number	Artwork Size H" x W" 1/4" bleed on all sides	Finished Size H" x W"
5DT-AOD3	5DT-32D3	12½" X 9¾"	12" X 91/4"

Custom Offer Sizes:

Gloss & Fabric Tech D		Tech Device Tattoos	Device Tattoos Matte Tech Device Tattoos	
Custom	Artwork Size H" x W" 1/4" bleed on all sides	Finished Size H" x W"	Artwork Size H" x W" 1/4" bleed on all sides	Finished Size H" x W"
Minimum	5½" x 5½"	5" X 5"	5½" x 5½"	5" X 5"
Maximum	59" x 198½"	58½" X 198"	53" x 198½"	52½" X 198"

^{*}For tech device tattoos under 8" x 8" minimum quantities may apply

Format/Layout Options **Available for Standard** Offer Sizes:



Artwork Requirements: Artwork size should include a 1/4" bleed on all sides. Bleed area should always be included in artwork even if design does not bleed.

CMYK processed file, Create full color process files with Adobe Acrobat PDF file, *note: if RGB or other non-CMYK colors are used

in your files the necessary conversions to CMYK will change the appearance of the colors when printed

Minimum 300 dpi, 600 dpi is recommended for highest quality output

File formats should be Tiffs, PDFs or Illustrator files with the fonts embedded in the file, *PDF files are preferred

Do not include any punctuation in the PDF file name and keep the name as short as possible

Order Process: Electronic or Manual

Proofs: Emailed Electronic PDF Proofs, add 2 working days to the item production. First proof is free of charge. Additional proofs

available for an additional charge. Production time begins after receipt of complete order including proof approval.

Production Samples, physical product proofs are available using customer artwork. These are offered at a % off the one

quantity price. Standard production times will apply.

Complimentary Proof may be required for certain orders to ensure satisfaction. Production time will be begin upon approval of

Packaging: Kraft tube or box shipping container depending on the size.

Multiple tech device tattoos will ship in one shipping container, multiple containers pending size and quantity.

Changes or Cancellations: Changes to orders in progress may require additional charges and/or time. Cancellations for orders in progress will require

payment for the order.